

You can use the Image Editor to display X-ray images and work with them. The Image Editor window is displayed automatically when you select to open an image. The Image Editor can be used to perform many operations on X-ray images. Using the editor, you can flip images vertically and horizontally, pan, zoom, change brightness and contrast settings, and much more.

Mithilfe der Bildbearbeitung können Sie Röntgenaufnahmen anzeigen und bearbeiten. Das Fenster für die Bildbearbeitung erscheint automatisch beim Öffnen eines Bildes. Die Bildbearbeitung bietet viele Möglichkeiten der Bearbeitung von Röntgenbildern. Sie können Röntgenbilder horizontal und vertikal spiegeln, schwenken, vergrößern, Helligkeit und Kontrast verändern und noch vieles mehr.